





Release Notes

Web 1.7.2
Mobile 2.2.7

Release Date:
11_2_22



New Features & Functions in this release

- Web (1.7.2) and Mobile (2.2.7)
 1. Changed compression to “Choose Any” from “Choose One” on Web App
 2. Added compression therapy question to Mobile App (“Choose Any”)
 3. Added offload pressure question to Mobile App (yes or no question)
 4. Added “Rotate Photo” feature (used when photo is taken with mobile device in a horizontal position) to the Mobile App

BUG FIXES and Optimizations

- none

1. Changed compression to “Choose Any” (aka multi-select) from “Choose One” on Web App

Last Name
Frederickson

First Name
Fred

Gender
Male

DOB
10/28/2011

Medical Record
-

✕

Wounds, Assessments & Treatments

Add / Edit Assessment

Wound
1

Assessment
1

Location
Torso Anterior (upper)

Assessment
(Required)

Debride, Cleanse / Wound
(Optional)

Wound Bed
(Optional)

Misc.
(Optional)

Risk Factors
(Optional)

Assessment Date
10/28/2022

Wound Measurements (centimeters)
Length 7.2 Width 7.1 Depth 7

Exudate/Drainage Amount (Choose One)

None

Light

Moderate

Heavy

Copious

Wound Thickness/Stage
Full Thickness

Wound Type
Venous Insufficiency Ulcer

Compression therapy

☐ None

☐ Inelastic Bandage

☐ Stockings or custom garment

☒ Leg elevation

☐ Elastic Bandage

☒ Dual Compression

☐ Pneumatic compression pump

*Describe any therapy detail in the assessment notes below (including # of layers)


Assessment Notes (optional)
Enter Notes

Save

Return

Treatment

2. Added compression therapy question to Mobile App (“Choose Any”)



Note: This is a required mobile app screen that only appears if the wound type = Venous Ulcer

3. Added offload pressure question to Mobile App (yes or no question)

The screenshot shows a mobile app interface for an assessment. At the top, the status bar displays the time 3:33, LTE signal, and battery level. Below the status bar is a green 'Back' button. The main heading is 'Answer assessment questions' with a clipboard icon and an 'Offline' status indicator. A progress bar shows six steps, with the fifth step (the current question) highlighted in blue. The question is: 'Device(s) being used to offload pressure and / or reduce sheer?'. Below the question are two input fields: 'Yes' and 'No'. At the bottom, there is a 'Note:' section with the text 'Describe device(s) used in the assessment notes' and a large text area for notes.

3:33 LTE

< Back

Answer assessment questions

Offline

Device(s) being used to offload pressure and / or reduce sheer?

Yes


No

Note:
Describe device(s) used in the assessment notes



Note: This is a required mobile app screen that only appears if the wound type = Pressure Ulcer

4. Added “Rotate Photo” feature (used when photo is taken with mobile device in a horizontal position) to the Mobile App


In this example (back heel of foot), photo was taken with patient laying on a horizontal surface. To properly take a wound in this situation you must turn the mobile device to a horizontal position to ensure the top of the phone is pointing at the head of the patient . .



Head of patient is →




Tap the “Rotate Photo” icon once to rotate the image 90 degrees




The photo will rotate and show on the screen like this . .

Continue as you would normally by measuring the reference marker and then the wound length and width . .



Remember: Length is measured “Head to Toe”



Remember: Width is measured “Hip to Hip”

